EAR SENSOR

Model CWS4

User Manual

English (UK)

English (US)

Nederlands

Deutsch

Française

Español (ES)

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
AHEAD OF THE HERD TOGETHER

CONTENTS

1	PRODUCT DESCRIPTION	3
1.1	Intended Use	3
1.2	Limitation of liability	3
1.3	Symbols and Markings	3
1.4	Main Parts	4
1.5	Typeplate	5
1.6	Technical Specifications	6

2	SAFETY WARNINGS	7
----------	------------------------	----------

3	INSTALLATION AND PLACEMENT	8
3.1	Installation Instructions	8
3.2	Preparing the Ear Sensor	9
3.3	Placing the Ear Sensor	10
3.4	Replacing the Ear Sensor	11

4	MAINTENANCE	11
----------	--------------------	-----------

5	TROUBLESHOOTING	12
----------	------------------------	-----------

6	DISPOSAL	13
----------	-----------------	-----------

1 PRODUCT DESCRIPTION

1.1 Intended Use

The CWS4 (further referred to as ear sensor) is part of the CowManager system and is meant to be installed in a cow's ear to collect information about the ear temperature, activity, rumination and eating behaviour of a cow. The device transmits the collected data via wireless communication to the system.

1.2 Limitation of liability

Installing, handling, and use of the product by the customer is beyond company control. Therefore, CowManager does not assume responsibility for loss, damage or expense resulting from improper installation, handling or misuse including – without limitation- bodily harm, injury, or damage to property, in connection with system information and the instructions set forth in this manual.

CowManager B.V. and its affiliates are not liable for damages and/or losses related to security breaches, any unauthorized access, interference, intrusion, leakage and/or theft of data or information, except for as explicitly provided by the terms of an agreement between the parties, see details in CowManager Terms and Conditions and NLDigital Terms.

1.3 Symbols and Markings

The approval overview can be found in the back of this booklet.

Disposal of the ear sensor; please do not throw away the ear sensor with the general household waste. Please use the designated separate collection system for those products (e.g. small electronics, small chemical waste).

IP67

Ingress protection code is 67.

The CE mark shows that the ear sensor meets the requirements for safety and health (Article 3.1(a)), electromagnetic compatibility (Article 3.1(b)), and the efficient use of the radio spectrum (Article 3.2) of the RED directive 2014/53/EU.

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.


Indicates useful information regarding the (use of the) product.

1.4 Main Parts

Accessories

- | | |
|-----------------|---------------------|
| 1 Insert ring | 5 Female tag |
| 2 LED indicator | 6 Male tag |
| 3 Ear sensor | 7 Activation magnet |
| 4 Tagger | |


CowManager part numbers

Ear Sensor	CM000R60 / CM000S90
Tagger	CM004003
Female tag	CM004012
Male tag	CM004011
Activation magnet	CM002004

1.5 Typeplate

Every ear sensor has a unique serial number. This number is located at the back of the ear sensor.

The ear sensor is available in two variants, the CWS4-S8 and CWS4-S9. Both ear sensors use the same electronics and have identical functionality. The only difference is the production method.

- ① Certification marks
 - ② Product details
 - ③ QR code
 - ④ Serial number
- i Remark:** The real serial number is longer than the two letters and six digits that are displayed on the back of the ear sensor. The number starts with "A0000ABB" but since that part is generic for all ear sensors, this has been omitted. If the QR-code is scanned all 16 positions of the serial number will be displayed.

1.6 Technical Specifications

Parameter	Value	Unit
Power		
Supply voltage	3.6	V
Battery capacity (at 20 °C / 68 °F)	1200	mAh
Battery chemistry	Lithium Thionyl	-
Wireless		
Standard	IEEE802.15.4a	-
Operational frequency	2405 / 2480	MHz
Operating Zigbee channel (default)	11	-
Modulation	OQPSK	-
Transmission power	<10	dBm
Antenna polarization	Omnidirectional	-
Range of RF link (outdoor, free line of sight)	100 / 330	m / ft
Environmental		
Environmental protection level / IP rating	IP67	-
Operating temperature	-20 – 40	°C
	-4 – 104	°F
Operating relative humidity	10 – 95	%RH
Storage temperature	0 – 40	°C
	32 – 104	°F
Physical		
Dimensions	50 × 65 × 19	mm
Weight	30	g

2 SAFETY WARNINGS

WARNING

- The ear sensor is meant to be used for the monitoring of cows only.
- Do not use the ear sensor on another position than the one indicated in **Chapter 3.3**.

Placing the Ear Sensor.


- Do not modify or alter the ear sensor in any way. Alteration or modification may affect the safety, performance and/or accuracy.
- Do not expose the ear sensor to fire or other heat sources.
- Avoid storage in direct sunlight.
- Do not crush or disassemble the ear sensor and do not subject the ear sensor to mechanical shock.
- To prevent damage, do not soak or immerse the ear sensor in any liquid solution.
- For cleaning of the ear sensor follow the instructions as described in **Chapter 4**.

Maintenance.

- Do not expose the ear sensor to excessive heat (above 80°C or 176°F). This may cause electrolyte leakage or explosion of the battery in the device. Don't expose the ear sensor to excessive cold (below -40°C or -40°F).

3 INSTALLATION AND PLACEMENT

3.1 Installation Instructions

To use a CowManager ear sensor, you must first activate the ear sensor in your CowManager system. For activation, the CowManager application must be installed on your PC, and a so-called “Coordinator” must be connected and configured.


Details regarding the installation of the CowManager system can be found in the CowManager System Installation Manual. Scan the QR code to see the complete CowManager System Installation Manual online.


It is important to correctly attach the ear sensor to the cow’s ear in a secure manner to avoid them from turning or falling off. The ear sensor can be attached to the ear either with the use of an ear sensor compatible EID or a blank tag.

Before using (electronic) ear tags for attaching the CowManager ear sensor, you should familiarize yourself with the local laws and regulations regarding the use of ear tags in cows.


3.2 Preparing the Ear Sensor

1. Put the new (female) ear tag (blank tag or RFID) in the ear sensor.


2. Remove the black insert of the tagger.

3. Place the ear sensor with the female tag and the male tag in the tagger.
-

3.3 Placing the Ear Sensor

i We advise only tagging calves older than 14 days and weighing more than 50 kg (110 lbs).

1. Disinfect the ear of the animal, before applying the ear tag and ear sensor to prevent infection. Use alcohol or chloric hexedine. Disinfection wipes (as used with dry-off) can be used as well.
 2. Place the ear sensor in the ear according to the illustration above.
- i** Make sure the ear sensor is always placed at the correct position in the ear. Choose the ear with the most space without other ear tags interfering, and do not attach the ear sensor at the same side as the moving part of the headlock.
3. Make sure that the correct ear sensor is linked to the correct animal. You can easily link the ear sensors to the animals by using the CowManager mobile app.

3.4 Replacing the Ear Sensor

Delink Ear Sensors:

1. Take the ear sensors off the ear of the cow, before marking them for replacement.
2. Delink the ear sensors in the CowManager application.

Replace:

Now that the ear sensors are delinked, mark them as non-functioning or lost.

1. Go to 'Input' > 'Replace Sensors'.
2. Click on the checkbox in front of the Sensor ID and choose 'Replace as a non-functioning sensor' or 'Replace as lost sensor'. A pop-up shows, tick the checkboxes and click 'Agree'.

NOTE! Only main users and users with access to 'Input' can replace ear sensors.

4 MAINTENANCE

The ear sensor is moulded with a waterproof material that also protects the ear sensor against barn gases and manure. If the ear sensor needs cleaning, proceed as follows:

- Always use water without additives, with a maximum temperature of 40°C (104°F), and a soft brush.
- Do not wash the ear sensor in a washing machine.
- The ear sensor should not be submerged.
- Do not dry the ear sensor in an oven or microwave but always use a cloth.
- Always handle ear sensors with care.

5 TROUBLESHOOTING

In order to operate correctly, the ear sensor needs to be able to communicate every 15 minutes with a CowManager router or coordinator. Communication can be disturbed because of the following reasons:

1. If the distance between the ear sensor and coordinator/router is too large. The operating distance of the ear sensor is around 100 meters (about 330 feet) in a free line of sight situation. Obstacles like walls, roofs, metal objects, foliage, geography or bad weather conditions will have a negative influence on the operating distance.
 -  Sensor data will be stored in the ear sensor for 14 days, before it will be overwritten.
2. A strong transmitter that operates in the same frequency band as the ear sensor. For example, Wi-Fi equipment like routers, access points and repeaters, but also P2P devices. The ear sensor communicates using Zigbee channel 11 in the 2,4GHz frequency band and this overlaps with channel 1 - 5 of a Wi-Fi network. In case of communication problems, check the configuration of Wi-Fi and other network equipment to make sure that this equipment uses Wi-Fi channel 6 or higher, minimizing or eliminating interference with CowManager devices.

When the environmental temperature is outside the range of $-20 - 40^{\circ}\text{C}$ / $-4 - 104^{\circ}\text{F}$ (See **Chapter 1.6. Technical Specifications**) the operating distance of the device may decrease and irreversible damage could occur. If the ear sensor cannot communicate because of a technical issue with the ear sensor itself, this is made visible by a blinking red LED. See **Chapter 1.4. Main Parts** for the location of the LED. Replace the ear sensor if this occurs.

6 DISPOSAL

The ear sensor is marked with this symbol. It means that electrical and electronic products should **not** be mixed with general household waste. There is a separate collection system for these products.

Attention: If you want to dispose of this ear sensor, please take this product to electronic waste collection points. Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling. Please contact your local authority for further details of your nearest electronic waste collection point.

CONTENTS

1	PRODUCT DESCRIPTION	15
1.1	Intended Use	15
1.2	Limitation of liability	15
1.3	Symbols and Markings	15
1.4	Main Parts	16
1.5	Typeplate	17
1.6	Technical Specifications	18
2	SAFETY WARNINGS	19
3	INSTALLATION AND PLACEMENT	20
3.1	Installation Instructions	20
3.2	Preparing the Ear Sensor	21
3.3	Placing the Ear Sensor	22
3.4	Replacing the Ear Sensor	23
4	MAINTENANCE	23
5	TROUBLESHOOTING	24
6	DISPOSAL	25
7	COMPLIANCE	25
7.1	FCC and IC Declarations	25

1 PRODUCT DESCRIPTION

1.1 Intended Use

The CWS4 (further referred to as ear sensor) is part of the CowManager system and is meant to be installed in a cow's ear to collect information about the ear temperature, activity, rumination and eating behavior of a cow. The device transmits the collected data via wireless communication to the system.

1.2 Limitation of liability

Installing, handling, and using the product by the customer is beyond the company's control. Therefore, CowManager does not assume responsibility for any loss, damage or expense resulting from improper installation, handling, or misuse including – without limitation – bodily harm, injury, or damage to property in connection with system information and the instructions set forth in this manual.

CowManager B.V. and its affiliates are not liable for damages and/or losses related to security breaches, any unauthorized access, interference, intrusion, leakage, and/or theft of data or information, except for as explicitly provided by the terms of an agreement between the parties, see details in CowManager Terms and Conditions and NLDigital Terms.

1.3 Symbols and Markings

The approval overview can be found in the back of this booklet.

Disposal of the ear sensor; please do not throw away the ear sensor with the general household waste. Please use the designated separate collection system for those products (e.g. small electronics, small chemical waste).

IP67

Ingress protection code is 67.

The CE mark shows that the ear sensor meets the requirements for safety and health (Article 3.1(a)), electromagnetic compatibility (Article 3.1(b)), and the efficient use of the radio spectrum (Article 3.2) of the RED directive 2014/53/EU.

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.


Indicates useful information regarding the (use of the) product.

1.4 Main Parts

Accessories

- | | |
|-----------------|---------------------|
| 1 Insert ring | 5 Female tag |
| 2 LED indicator | 6 Male tag |
| 3 Ear sensor | 7 Activation magnet |
| 4 Tagger | |


CowManager part numbers

Ear Sensor	CM000R60 / CM000S90
Tagger	CM004003
Female tag	CM004012
Male tag	CM004011
Activation magnet	CM002004

1.5 Typeplate

Every ear sensor has a unique serial number. This number is located at the back of the ear sensor.

The ear sensor is available in two variants, the CWS4-S8 and CWS4-S9. Both ear sensors use the same electronics and have identical functionality. The only difference is the production method.

- ① Certification marks
 - ② Product details
 - ③ QR code
 - ④ Serial number
- i Remark:** The real serial number is longer than the two letters and six digits that are displayed on the back of the ear sensor. The number starts with "A0000ABB" but since that part is generic for all ear sensors, this has been omitted. If the QR-code is scanned all 16 positions of the serial number will be displayed.

1.6 Technical Specifications

Parameter	Value	Unit
Power		
Supply voltage	3.6	V
Battery capacity (at 20 °C / 68 °F)	1200	mAh
Battery chemistry	Lithium Thionyl	-
Wireless		
Standard	IEEE802.15.4a	-
Operational frequency	2405 / 2480	MHz
Operating Zigbee channel (default)	11	-
Modulation	OQPSK	-
Transmission power	<10	dBm
Antenna polarization	Omnidirectional	-
Range of RF link (outdoor, free line of sight)	100 / 330	m / ft
Environmental		
Environmental protection level / IP rating	IP67	-
Operating temperature	-20 – 40	°C
	-4 – 104	°F
Operating relative humidity	10 – 95	%RH
Storage temperature	0 – 40	°C
	32 – 104	°F
Physical		
Dimensions	50 × 65 × 19	mm
Weight	30	g

2 SAFETY WARNINGS

WARNING

- The ear sensor is meant to be used for the monitoring of cows only.
- Do not use the ear sensor on another position than the one indicated in **Chapter 3.3**.

Placing the Ear Sensor.


- Do not modify or alter the ear sensor in any way. Alteration or modification may affect the safety, performance and/or accuracy.
- Do not expose the ear sensor to fire or other heat sources.
- Avoid storage in direct sunlight.
- Do not crush or disassemble the ear sensor and do not subject the ear sensor to mechanical shock.
- To prevent damage, do not soak or immerse the ear sensor in any liquid solution.
- For cleaning of the ear sensor follow the instructions as described in **Chapter 4**.

Maintenance.

- Do not expose the ear sensor to excessive heat (above 80°C or 176°F). This may cause electrolyte leakage or explosion of the battery in the device. Don't expose the ear sensor to excessive cold (below -40°C or -40°F).

3 INSTALLATION AND PLACEMENT

3.1 Installation Instructions

To use a CowManager ear sensor, you must first activate the ear sensor in your CowManager system. For activation, the CowManager application must be installed on your PC, and a so-called “Coordinator” must be connected and configured.


Details regarding the installation of the CowManager system can be found in the CowManager System Installation Manual. Scan the QR code to see the complete CowManager System Installation Manual online.


It is important to correctly attach the ear sensor to the cow’s ear in a secure manner to avoid them from turning or falling off. The ear sensor can be attached to the ear either with the use of an ear sensor compatible EID or a blank tag.

Before using (electronic) ear tags for attaching the CowManager ear sensor, you should familiarize yourself with the local laws and regulations regarding the use of ear tags in cows.


3.2 Preparing the Ear Sensor

1. Put the new (female) ear tag (blank tag or RFID) in the ear sensor.


2. Remove the black insert of the tagger.

3. Place the ear sensor with the female tag and the male tag in the tagger.
-

3.3 Placing the Ear Sensor

i We advise only tagging calves older than 14 days and weighing more than 50 kg (110 lbs).

1. Disinfect the ear of the animal before applying the ear tag and ear sensor in order to prevent infection. Use alcohol or chloric hexedine. Disinfection wipes (as used with dry-off) can be used as well.
 2. Place the ear sensor in the ear according to the illustration above.
- i** Make sure the ear sensor is always positioned correctly in the ear. Choose the ear with the most space, without other ear tags interfering, and do not attach the ear sensor on the same side as the moving part of the headlock.
3. Ensure the correct ear sensor is linked to the correct animal. You can easily link the ear sensors to the animals by using the CowManager mobile app.
-

3.4 Replacing the Ear Sensor

Delink Ear Sensors:

1. Remove the ear sensors from the cow's ear before marking them for replacement.
2. Delink the ear sensors in the CowManager application.

Replace:

Now the ear sensors are delinked, mark them as non-functioning or lost.

1. Go to 'Input' > 'Replace Sensors'.
2. Click on the checkbox in front of the Sensor ID and choose 'Replace as a non-functioning sensor' or 'Replace as lost sensor'. A pop-up appears, tick the checkboxes and click 'Agree'.

NOTE! Only main users and users with access to 'Input' can replace ear sensors.

4 MAINTENANCE

The ear sensor is moulded with a waterproof material that also protects the ear sensor from barn gases and manure. If the ear sensor needs cleaning, proceed as follows:

- Always use water without additives, with a maximum temperature of 40°C (104°F), and a soft brush.
- Do not wash the ear sensor in a washing machine.
- The ear sensor should not be submerged.
- Do not dry the ear sensor in an oven or microwave but always use a cloth.
- Always handle ear sensors with care.

5 TROUBLESHOOTING

In order to operate correctly, the ear sensor needs to be able to communicate every 15 minutes with a CowManager router or coordinator. Communication can be disturbed because of the following reasons:

1. If the distance between the ear sensor and the coordinator/router is too large. The operating distance of the ear sensor is around 100 meters (approximately 330 feet) in a free line of sight situation. Obstacles such as walls, roofs, metal objects, foliage, geography, or adverse weather conditions will negatively affect the operating distance.
 -  Sensor data will be stored in the ear sensor for 14 days, before it will be overwritten.
2. A strong transmitter that operates in the same frequency band as the ear sensor. For example, Wi-Fi equipment like routers, access points and repeaters, but also P2P devices. The ear sensor communicates using Zigbee channel 11 in the 2.4GHz frequency band and this overlaps with channel 1 - 5 of a Wi-Fi network. In case of communication problems, check the configuration of Wi-Fi and other network equipment to make sure that this equipment uses Wi-Fi channel 6 or higher, minimizing or eliminating interference with CowManager devices.

When the environmental temperature is outside the range of $-20 - 40^{\circ}\text{C}$ / $-4 - 104^{\circ}\text{F}$ (See **Chapter 1.6. Technical Specifications**) the operating distance of the device may decrease and irreversible damage could occur. If the ear sensor cannot communicate because of a technical issue with the ear sensor itself, this is made visible by a blinking red LED. See **Chapter 1.4. Main Parts** for the location of the LED. Replace the ear sensor if this occurs.

6 DISPOSAL

The ear sensor is marked with this symbol. It means that electrical and electronic products should **not** be mixed with general household waste. There is a separate collection system for these products.

Attention: If you want to dispose of this ear sensor, please take this product to electronic waste collection points. Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling. Please contact your local authority for further details of your nearest electronic waste collection point.

7 COMPLIANCE

7.1 FCC and IC Declarations

Compliance statement (part 15 . 19)

This device complies with part 15 of the FCC Rules and with ISDE Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions:

1. this device may not cause harmful interference, and
2. this device must accept any interference received, including interference that may cause undesired operation.

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Warning (part 15. 21)

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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Approvals

Certification approval logo can be found on the product (see **Chapter 1.3. Symbols and Markings**). In case of limited space, certification approval logo can be found on packaging according to regulations.

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